

The refined choice



The Fine Range - accuracy you demand and quality as standard

Complimenting the traditional Swann-Morton blade shapes and developed potentially for more complex and challenging surgical procedures which may be encountered in such disciplines as reconstruction, ophthalmics, cardiology and even hair restoration, the surgeons may sometimes prefer our Fine Range.

The Fine Range has small, slim, precise blades in stainless steel for delicate operations.

Six special handles, including the SF13 and SF23 with independently screwed top sections to

facilitate blade-changing, take all blades, in stainless steel, sterile, individually sealed in metal foil packets and boxed in units of 25.

Handles are supplied individually in plastic wallets.

Sterile blades should be selected for all invasive surgical procedures and are for single patient use thus minimising the risk of cross infection.

Swann-Morton®
Made in Sheffield - *Precisely*



CE
All medical devices are CE marked in accordance with the Medical Device Directive (93/42/EEC)



Swann-Morton®

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

Fine Range Blades - Sterile

Blades in Stainless Steel - Sterile, individually V.P.I. protected in Metal Foil Packets

Fits handles SF1, SF2, SF3, SF4, SF13 & SF23

| Product | Product No | Unit | Product | Product No | Unit |
|---------|------------|--------------|-------------|------------|-------------------------------------|
| SM61 | 5901 | | SM61(SB) | 5911 | |
| SM62 | 5902 | | SM62(SB) | 5912 | |
| SM64 | 5904 | 25 blades | SM63 | 5903 | 25 blades |
| SM65 | 5905 | individually | SM69 | 5909 | individually |
| SM65A | 5906 | packed | SP90 | 5921 | packed |
| SM67 | 5907 | | SP91 | 5922 | |
| SM68 | 5908 | | Myringotomy | 2102 | 10 blades individually packed |

Specialist Applications

Cardio-Vascular Surgery: SM65, SM65a, SM67. Plastic Surgery: SM63, SM64, SM67 & SM69. Ophthalmology: SM63, SM64, SM67 & SM69. Orthopaedic Surgery: SM64, SM65 & SM65a. E.N.T. Surgery: SM68 & Myringotomy. Podiatry: SM61, SM61(SB), SM62, SM62(SB), PD81 & PD82. Dental Periodontics & Implantology: SM65, SM67 & SM69. Hair Restoration & Transplantation: SP90 & SP91

Fine Range Handles - Non-Sterile

Handles in Stainless Steel - Non-Sterile, individually sleeved

1 handle individually packed (2 piece)

| Product | Product No | Unit |
|---------|------------|------|
| SF1 | 6051 | 1 |
| SF2 | 6052 | 1 |
| SF3 | 6053 | 1 |
| SF4 | 6054 | 1 |

1 handle individually packed (3 piece)

| Product | Product No | Unit |
|---------|------------|------|
| SF13 | 6061 | 1 |
| SF23 | 6062 | 1 |



Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.